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
## Set up a sequencing run

1. Navigate to the Start page and click **Start Sequencing**.

2. In the Positions tab, fill in the experiment name, sample ID, and select the flow cell type. Click **Continue to kit selection**.


### 2 Type in the experiment name, sample ID and choose flow cell type from the drop down menu.

Ensure the experiment name and sample ID does not contain any personally identifiable information.

**Note:** If sample ID is not filled in, there will be no sample ID in the folder structure.

To use previously saved settings, please see the "Save sequencing run settings" page of this protocol.

Select **Continue** to move to the Kit Selection page.

3. Select the kit and any expansion(s) used to prepare the library. Click **Continue to run options**.

### 3 Select the kit and expansion(s) used to prepare the library.

One sequencing kit can be selected as well as one or more expansions. If a sequencing kit does not have any associated expansions, no expansion options will appear.

The filter options may be used to find the kit used. For example, click **DNA** to filter for DNA sequencing kits.

Click **Continue** to navigate to the run configuration page.

Sample type	PCR-free	Multiplexing	Control
<input type="radio"/> DNA	<input type="radio"/> PCR	<input type="radio"/> Yes	<input type="radio"/> Control
<input type="radio"/> RNA	<input type="radio"/> PCR-free	<input type="radio"/> No	<input type="radio"/> Control

Reset filters

Sequencing Kit	Kit ID	Selected
Ligation Sequencing Kit	SOK-LSK114	<input type="radio"/>
Native Barcoding Sequencing Kit 96	SOK-NB0114-96	<input checked="" type="radio"/>
PCR cDNA Barcoding Kit - 24	SOK-PCB114-24	<input type="radio"/>
Rapid PCR Barcoding Kit	SOK-RPB114-24	<input type="radio"/>
Ligation Sequencing Kit	SOK-LSK114-3L	<input type="radio"/>
Ultra-Long DNA Sequencing Kit	SOK-ULK114	<input type="radio"/>
16S Barcoding Kit	SOK-16S114-24	<input type="radio"/>
Rapid Sequencing Kit	SOK-RAD114	<input type="radio"/>
Native Barcoding Sequencing Kit 24	SOK-NB0114-24	<input checked="" type="radio"/>
Multiplex Ligation Sequencing Kit XL	SOK-MLK114-96-3L	<input type="radio"/>
Rapid Barcoding Kit 24	SOK-RBK114-24	<input type="radio"/>

Select expansion pack

Expansion Pack	Expansion ID	Selected
Cas9 Barcoding Expansion	EXP-CB0001	<input checked="" type="radio"/>

Back Continue

4. Specify the sequencing run length and minimum read length or keep the default settings. Click **Continue to analysis**.

5. Choose your basecaller model and select any available barcoding and alignment options or keep the default settings. Click **Continue to output**.


### 5 Choose your basecalling options:

The basecalling option is on by default to perform real-time basecalling. This can be switched off and basecalling performed post-sequencing.

**Note:** If basecalling is disabled, barcoding and alignment will not be performed during sequencing.

Modified bases is off by default. This can be switched on to use the CpG context models and basecall 5mC and 5hmC.

To specify your basecall model and modified bases, click the **Edit** icon to open basecalling options dialogue box.

## Basecall models

Basecalling models can be selected at two stages in the MinKNOW software:

- **Real-time basecalling:** A basecalling model can be selected when setting up a sequencing run on the "Run configuration" page of MinKNOW. This enables basecalling of the sequencing run in real-time.
- **Post-run basecalling:** Basecalling can be set up post-sequencing run on MinKNOW. Instructions can be found in the "[Post-run basecalling](#)" section of this protocol.

Dorado is the basecaller server in MinKNOW which provides multiple models for basecalling nanopore data.

There are three options for model selection in the drop down menu:

**1. High-accuracy basecalling (HAC)** - This is our default recommended basecalling model that provides a high single-molecule accuracy. It is currently 5-8 times more computationally-intensive than the Fast model, so users should ensure their data transfer and device utilisation is scaled appropriately for this.

**2. Super-accurate basecalling (SUP)** - The Super accurate model has an even higher single-molecule accuracy, and is ~3 times more intensive than the HAC model.

**3. Fast basecalling** - This model is able to keep up with a high-throughput sequencing experiment on a MinION, GridION or PromethION.

Note: GPU devices will basecall in near real-time; CPU devices, like a standard laptop for MinION Mk1B/Mk1D, will not maintain real-time basecalling.

6. Choose your barcoding option

**6 Choose your barcoding options:**

Barcoding options are only available when a barcoding sequencing kit or expansion has been selected. Barcoding can be switched off and performed post-sequencing.


To specify the barcoding options, click the **Edit** icon to open the barcoding options dialogue box. From here, barcode trimming and custom barcode selection can be defined.

By default, barcoding trimming is off. This can be enabled to remove barcodes from demultiplexed reads. Please note, some primer sequences may also be trimmed together with the barcodes.

By default, custom barcode selection is off. This can be enabled to only select for specified barcodes used in your library preparation for the current sequencing run.

7.

Specify your output data location, format and filtering options or keep the default settings. Click **Continue to final review**.

8. Click **Start**.